




BOT 3000 Measurement Results



<b>Method:</b>	Dynamic Coefficient of Friction	<b>Contact Time:</b>	-	<b>Distance:</b>	1 ft				
<b>Device:</b>	BOT-3000 1053 SlipDoctors	<b>Calibration Date:</b>	Wed 07/22/09	<b>Verification Date:</b>	Thu 10/01/09				
<b>Slider:</b>	10267 <b>Type:</b> Neolite®	<b>Mfg. Date:</b>	Tue 04/14/09	<b>No. of Measures:</b>	113	<b>Total Distance:</b>	67 ft		
<b>Result:</b>	<b>Date:</b> Thu 10/08/09 02:10 PM	<b>Min. Value:</b>	0.33	<b>Avg. Value:</b>	0.35	<b>Max. Value:</b>	0.40	<b>Facility Code:</b>	----